F-834 User Manual


[Function]

- 1. Eliminate the aging stratum corneum, tender skin, clear away dirt from pores, fade acne scars;
- 2. Suck accumulated oil from the skin and the pore, improve dull skin, eliminate facial swelling, eye swelling;
- 3. Increase skin elasticity, eliminate slender wrinkles, improve skin texture, keep the skin youth;
- 4. Improve the skin from dull, rough, uneven complexion, pigment, spots;
- 5. Dilute the tattoo pigment, fat pattern, stretch marks;
- 6. Eliminate the aging stratum corneum, tender skin.

[Fundamental]

This Instrument is made use of a mosaic of different thickness diamond particles, with controlled grating, and matched up the vacuum suction to control strength. It makes partial miniature rejuvenation for the lymphatic drainage, skin aging, pigmentation or large pores, etc and brings the most gentle and rapid regeneration for the skin. Not only can accelerate the renewal of the aging stratum corneum, but also increase the full absorption of skin care products.

[manipulate]



- 1 power light 2 mains switch
- 3 Suction regulator 4 Negative pressure indicator

[Features]

- 1. The newest high-tech diamond dermabrasion rejuvenation technology, security, painless, no side-effects:
- 2. 9, Diamond Head: 3 pins for different facial skin, 2 specialty pins for eye and slender wrinkles, 4, specialty pins for the body with a wide range professional treatment
- 3. Humanization design, developed operating handle elaborately, accord with hand structure and facilitate the beauticians;
- 4. With simplicity of operation, it makes the skin silky, natural, safe effectively;
- 5. No any powder in the process, to avoid the powder residue on the skin or suck by the guests or beautician;
- 6. Needn't rest after treatment, you can go out or make-up, it wouldn't affect the normal social life.